



Pipetting Proficiency Test Method




Background

- Pipetting variance can lead to errors which can:
 - Result in unrecognized wrong results
 - Failure to pass a good product (as in QC)
 - Failure to fail a poor product (as in QC)
- Pipetting variance can be attributed to:
 - Poorly seated tips
 - Variable viscosities not taken into account
 - Profusion (residual volume on exterior of tip) carry over
 - Pipette out of calibration (dropped or overly used)
 - Inexperience or need for retraining



Definition


Despite training or prior knowledge, it is possible to drift from uniform practice (as in a golf swing). While pipetting may be effective, it may not yield an accurate and/or consistent result.





Purpose

- This Test Method is designed to characterize pipetting efficiency across the organization.
- Depending on results, there may or may not be a need for training/re-training.
- Uniform pipetting practices are essential to the success of Salimetrics and to our customers.
- This test method, or a revision, is intended to be developed for our customers to help them ensure that their pipetting methods are correct.



Approach

- Who: All individuals who routinely operate pipettes should consider participating in this evaluation.
- What: Fluids of varying viscosity should be pipetted by each operator into microwell plates.
- When: Before performing assays. The procedure should only take a half hour.
- Where: At your own lab bench, with your current pipette
- Why – Ibidem. Make sure your results are consistent.



Methods

- Fluids of varying viscosities (xanthan gum, water, EtOH) containing yellow dye were developed for this study
- The fluids will be pipetted using different volumes in multiple replicates by each operator into microwell plates
- All additions are first pass only – no repeats

Single channel pipetting

10 uL (viscous) into dry wells, add 100 uL water, 48 replicates

50 uL (viscous) into dry wells, add 50 uL water, 48 replicates

100 uL (viscous) into dry wells, 48 replicates

Multi channel pipetting

50 uL (viscous) into dry wells, add 50 uL water, 48 replicates

100 uL (viscous) into dry wells, 48 replicates

200 uL (25% EtOH) into dry wells, 48 replicates



Materials Needed

Each operator should be provided with:

4 dye fluids made with McCormick's Yellow Dye

- 10 uL fluid (1 mL), dilute dye 58-fold in xanthan gum solution
- 50 uL fluid (8 mL), dilute dye 291-fold in xanthan gum solution
- 100 uL fluid (15 mL), dilute dye 580-fold in xanthan gum solution
- 200 uL fluid (15 mL), dilute dye 1160-fold in 25% EtOH solution

6 uncoated microwell plates

Xanthan Gum Solution

0.95 g/L Xanthan gum in water


0.35 g/L methyl p-hydroxy benzoate (preservative)

25% EtOH Solution

25% EtOH in water

0.35 g/L methyl p-hydroxy benzoate (preservative)

Pipetting Patterns – Single channel



10sc (single channel)

Add 10 uL dye fluid to each of well dry wells.
Add 100 uL water to each column of wells with
multichannel pipetter.
Read absorbance at 450 – 492.

	1	2	3	4	5	6
A						
B						
C						
D						
E						
F						
G						
H						

ADD
10 uL
Dye Fluid
→
to each
well

	1	2	3	4	5	6
A	↓	↓	↓	↓	↓	↓
B	↓	↓	↓	↓	↓	↓
C	↓	↓	↓	↓	↓	↓
D	↓	↓	↓	↓	↓	↓
E	↓	↓	↓	↓	↓	↓
F	↓	↓	↓	↓	↓	↓
G	↓	↓	↓	↓	↓	↓
H	↓	↓	↓	↓	↓	↓

single channel addition

ADD
100 uL
water
→
to each
well

	1	2	3	4	5	6
A	↑	↑	↑	↑	↑	↑
B						
C						
D						
E						
F						
G						
H						

multi channel addition

50sc (single channel)

Add 50 uL dye fluid to each of well dry wells.
Add 50 uL water to each column of wells with multichannel
pipetter. Read absorbance at 450 – 492.

	1	2	3	4	5	6
A						
B						
C						
D						
E						
F						
G						
H						

ADD
50 uL
Dye Fluid
→
to each
well

	1	2	3	4	5	6
A	↓	↓	↓	↓	↓	↓
B	↓	↓	↓	↓	↓	↓
C	↓	↓	↓	↓	↓	↓
D	↓	↓	↓	↓	↓	↓
E	↓	↓	↓	↓	↓	↓
F	↓	↓	↓	↓	↓	↓
G	↓	↓	↓	↓	↓	↓
H	↓	↓	↓	↓	↓	↓

single channel addition

ADD
50 uL
water
→
to each
well

	1	2	3	4	5	6
A	↑	↑	↑	↑	↑	↑
B						
C						
D						
E						
F						
G						
H						

multi channel addition

100sc (single channel)


Add 100 uL dye fluid to each of well dry wells.
Read absorbance at 450 – 492.

	1	2	3	4	5	6
A						
B						
C						
D						
E						
F						
G						
H						

ADD
100 uL
Dye Fluid
→
to each
well

	1	2	3	4	5	6
A	↓	↓	↓	↓	↓	↓
B	↓	↓	↓	↓	↓	↓
C	↓	↓	↓	↓	↓	↓
D	↓	↓	↓	↓	↓	↓
E	↓	↓	↓	↓	↓	↓
F	↓	↓	↓	↓	↓	↓
G	↓	↓	↓	↓	↓	↓
H	↓	↓	↓	↓	↓	↓


single channel addition



Pipetting Patterns - Multichannel


50mc (multi channel)

Add 50 uL dye fluid to each of well dry wells with all eight channels, left to right. Add 50 uL water to each column of wells with multichannel pipetter. Read absorbance at 450 – 492.




100mc (multi channel)

Add 100 uL dye fluid to each of well dry wells with all eight channels, left to right. Read absorbance at 450 – 492.



200mc (multi channel)

Add 200 uL dye fluid to each of well dry wells with all eight channels, left to right. Read absorbance at 450 – 492.





Analysis (GAGE R&R)

Uploaded results into Excel or equivalent analysis software program and assess for ANOVA

- well-to-well (repeatability)
- operator-to-operator (reproducibility)

**Note: Manual and automatic pipetter results may need to be handled separately.*



Example Evaluation


Methods

- Artificial saliva fluids with yellow dye were developed to deliver an OD of 1.0 after pipetting 10, 50, 100 or 200 uL into microwell plates.
- Plates were read on a plate reader and the data uploaded to a statistics program.
- Data was analyzed, by MiniTab: ANOVA, Box/Whisker, Dotplot and GAGE R&R

Example – single channel

10sc (single channel)

Add 10 uL dye fluid to each of well dry wells.
Add 100 uL water to each column of wells with multichannel pipetter.
Read absorbance at 450 – 492.






Example Evaluation

Example – multi channel

100mc (multi channel)

Add 100 uL dye fluid to each of well dry wells with all eight channels, left to right.

Read absorbance at 450 – 492.



Results

Thirteen (13) operators each pipetted dye solutions into six plates (48 wells each) for a total of 288 wells

- 10 uL, 50 uL, 100 uL single channel
- 50 uL, 100 uL, 200 uL multichannel

10 uL, single channel, worst case fluid, viscous, high dye concentration, small volume






<u>Within Operator</u>	
Avg CV =	8.2%
Stdev =	2.6%

<u>Operator - Operator</u>	
CV =	8%

<u>Pooled (Total)</u>	
CV =	11%

<u>Nested ANOVA:</u> 10 uL SC –vs. Operator, Replicate Variance Components	
Source	
Operator-Operator	46 %
Well-Well	54 %



50 uL, single channel, viscous

Within Operator

Avg CV = 8.3%

Stdev = 2.0%

Operator - Operator

CV = 9%

Pooled (Total)

CV = 12%



Nested ANOVA:

10 uL SC –vs. Operator, Replicate
Variance Components

Source

Operator-Operator 53 %

Well-Well 47 %



100 uL, viscous, single channel

Within Operator

Avg CV = 2.0%

Stdev = 1.0%

Operator - Operator

CV = 4%

Pooled (Total)

CV = 4%



Nested ANOVA:

10 uL SC –vs. Operator, Replicate
Variance Components

Source

Operator-Operator 75 %

Well-Well 25 %



50 uL, viscous, multi channel

Within Operator

Avg CV = 7.6%
Stdev = 2.2%

Operator - Operator

CV = 9%

Pooled (Total)



CV = 12%

Nested ANOVA:

10 uL SC –vs. Operator, Replicate
Variance Components

Source

Operator-Operator 58 %
Well-Well 42 %



100 uL, viscous, multi channel

Within Operator

Avg CV = 1.7%
Stdev = 0.7%

Operator - Operator

CV = 1.7%

Pooled (Total)



CV = 2.5%

Nested ANOVA:

10 uL SC –vs. Operator, Replicate
Variance Components

Source

Operator-Operator 44 %
Well-Well 56 %



200 uL, low viscosity, multi channel



Within Operator

Avg CV = 1.1%
Stdev = 0.3%

Operator - Operator

CV = 0.5%

Pooled (Total)

CV = 1.3%


Nested ANOVA:

10 uL SC –vs. Operator, Replicate
Variance Components


Source

Operator-Operator 11 %
Well-Well 89 %

200 uL Multi Channel



200 uL Multi Channel



Summary and Conclusions

- Levels of error are within expected ranges for manual pipetting.
- Larger volumes typically yield better precision. This is true in automated systems also.
- Pipette types (Hamilton, Finnpipette, Finnpipette 2, BioHit) were not factors (data not shown)
- Testing method was easy to implement, follow and execute.
- Remind operators to ensure tips are well seated and to minimize profusion (already common practice).